DE	Verarbeitungsanweisung für uvex medicare Einlagen und Komponenten zur Anfertigung
EN	Processing instructions for uvex medicare inserts and components for manufacturing
FR	Instructions d'utilisation pour les semelles uvex medicare et les composants de fabrication
IT	Istruzioni di lavorazione per solette uvex medicare e componenti di produzione
ES	Instrucciones de trabajo para plantillas ortopédicas uvex medicare y componentes para su fabricación
PT	Instruções de tratamento de palmilhas e componentes de fabrico uvex medicare
NL	Verwerkingsinstructie voor uvex medicare inlegzolen en componenten voor de vervaardiging
LT	Apdorojimo instrukcijos „uvex medicare“ vidpadžiams ir gamybos komponentams
LV	Apstrādes norādījumi uvex medicare ieliktņiem un izgatavošanas komponentiem
DK	Forarbejdningsanvisning til uvex medicare-indlæg og fremstillingskomponenter
NO	Tilpassingsinstruksjoner for uvex medicare-innlegg og -komponenter for tilpassing
EE	Töötlemisjuhised uvex medicare sisetaldadele ja komponendid nende valmistamiseks
CZ	Návod na zpracování pro vložky uvex medicare a komponenty pro výrobu
PL	Instrukcja dotycząca obróbki wkładek uvex medicare i komponentów do wykonania
SE	Bearbetningsanvisning för uvex medicare-iläggssulor och komponenter för tillverkning
FI	Käsittelyohje uvex medicare -tukupohjallisille ja komponenteille valmistamista varten
BG	Инструкции за обработка за стелки и производствени компоненти uvex medicare
HR	Upute o načinu obrade za uloške uvex medicare i komponente za izradu
RS	Uputstvo za obradu za uvex medicare uloške i komponente za proizvodnju
CN	处理说明 用于制造 uvex 医疗护理鞋垫和组件
HU	Kidolgozási utasítások az uvex medicare talpbetétek és alkatrészek gyártásához
UA	Інструкція з обробки устілок і компонентів uvex medicare для виготовлення
HI	प्रोसेस से जुड़े निर्देश uvex medicare के ईसर्ट्स और मैनुफैक्चरिंग से जुड़े कॉम्पोनेंट्स के लिए
JP	加工手順 製造向け uvex medicare の インサートとコンポーネント
SI	Navodila za obdelavo vložkov uvex medicare in komponent za izdelavo
RO	Instrucțiuni de prelucrare pentru inserții uvex medicare și componente pentru executare
RU	Инструкция по обработке супинаторов uvex medicare и компонентов для их изготовления
SK	Pokyny na spracovanie pre vložky uvex medicare a komponenty na výrobu
GR	Οδηγίες επεξεργασίας για υλικά και ενθέματα uvex medicare προς παραγωγή
TR	uvex medicare tabanlıklar ve üretim bileşenleri için işleme talimatları
AR	تعليمات المعالجة: لعمليات التصنيع الخاصة بالحشوات والمواد ذات الصلة بمنتجات uvex medicare

Weitere Information

More information

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www.uvex-safety-shop.de/home

Processing instructions

for uvex medicare inserts and
components for manufacturing

All components for manufacturing uvex comfortable customised inserts are tested and certified in combination with uvex safety footwear* according to EN ISO 20345. As regards product liability, please note that the type examination is only valid if materials approved by uvex are used and the processing instructions are followed.

You then act as the manufacturer towards your customer and are authorised to issue a declaration of conformity in your name. Declaration of conformity templates are available to download at:

<https://www.uvex-safety.com/en/knowledge/certificates/eu-declarations-of-conformity/>

uvex comfortable customised inserts consist of an EVA support core, a cover made from conductive PU and a perforated microfibre layer. The following materials can be used to make specific adjustments according to the indication:

Padding material: black, thickness: 7 mm (article number: 6497007)

Lifting material: black EVA, thickness: 8 mm (article number: 6497008)

In addition to the certified materials, any reinforcement material and/or pads can be used in line with the indication.

The cover is available separately as a replacement and for correcting customised inserts:

Cover approx. 350 x 450 mm + 2 pieces. uvex medicare logo (item number 6497009)

Please use the outer shoe material only in combination with the EVA support core for uvex medicare inserts.

Machining

Adjustments can be made using thermoplastic deformation at a temperature of approximately 80°C (warmed through). When heating specific points with hot air, do not exceed the three-second exposure time at an airflow temperature of 300°C. Reduce the exposure time accordingly when working at airflow temperatures above 300°C. The bottom length of the support core (preassembled insert and two-part insert set) depends on the insole and ends 1 cm before the protective toe cap starts. This is very important for maintaining the respective protective function of the shoe. Please pay particular attention when milling the support core out of the milling block (6497401) yourself. Changing the insert (of the top cover) in the area of the protective toe cap is not permitted and there must be no material residue, e.g. adhesive. A heel lift (max. 7 mm) is permissible in the heel area. The lateral height of the supporting mould can be freely selected.

Bonding

Commercially available polychloroprene or polyurethane adhesives can be used to bond the components. Make sure that all parts are free of dust and that there is no residue in the area of the protective toe cap. Put the components together and press at a pressure of 2 bar for approximately three seconds. When attaching a forefoot pad, make sure that the pad ends in the area of the protective toe cap and is bonded with antistatic adhesive.

Sanding

There must be no EVA material fitted/remaining in the area of the protective toe cap. The insert can be sanded afterwards — the following sanding belt grains are recommended:

Rough sanding: 24–40

Fine sanding: 80–120

* all uvex safety footwear suitable for uvex orthopaedic comfortable customised inserts is identified with the uvex medicare symbol in the uvex foot protection catalogue.

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